


Tracking


Thierry Pécot
Research Engineer
Biosit SFR UMS CNRS 3480 - Inserm 018
CZI Imaging Scientist




Tracking main steps



Frame #1




Frame #2




Frame #3

Tracking main steps


Segmentation



Frame #1




Frame #2




Frame #3

Tracking main steps

Data association, particle linking,...



Tracking main steps



TrackMate



Methods

Volume 115, 15 February 2017, Pages 80–90



TrackMate: An open and extensible platform for single-particle tracking

Jean-Yves Tinevez ^{a,*,1}, Nick Perry ^{a,1}, Johannes Schindelin ^{b,2}, Genevieve M. Hoopes ^c, Gregory D. Reynolds ^c, Emmanuel Lapointe ^d, Sebastian Y. Bednarek ^c, Spencer L. Shorte ^a, Kevin W. Eliceiri ^{b,*}

Bringing TrackMate into the era of machine-learning and deep-learning

^a Dmitry Ershov, Minh-Son Phan, ^b Joanna W. Pylvänäinen, ^c Stéphane U. Rigaud, Laure Le Blanc, ^d Arthur Charles-Orszag, ^e James R. W. Conway, ^f Romain F. Laine, ^g Nathan H. Roy, ^h Daria Bonazzi, ⁱ Guillaume Duménil, ^j Guillaume Jacquemet, ^k Jean-Yves Tinevez

doi: <https://doi.org/10.1016/2021.09.03.458852>

This article is a preprint and has not been certified by peer review [what does this mean?].

nature methods

Explore content ▾ About the journal ▾ Publish with us ▾ Subscribe

[nature](#) > [nature methods](#) > [articles](#) > [article](#)

Published: 20 July 2008

Robust single-particle tracking in live-cell time-lapse sequences

[Khuloud Jagaman](#) , [Dinah Loerke](#), [Marcel Mettlen](#), [Hirokata Kuwata](#), [Sergio Grinstein](#), [Sandra L. Schmid](#) & [Gaudenz Danuser](#)

TrackMate

1 - SEGMENTATION

Classic detectors:


- Manual annotation
- Threshold detector
- DoG/LoG detector



TrackMate

1 - SEGMENTATION


LoG detector: scale-based detector, powerful with objects of the same size such as intracellular particles



TrackMate

1 - SEGMENTATION



LoG detector: scale-based detector, powerful with objects of the same size such as intracellular particles



TrackMate

1 - SEGMENTATION

LoG detector: scale-based detector, powerful with objects of the same size such as intracellular particles



TrackMate

1 - SEGMENTATION

Classic detectors:

- Manual annotation
- Threshold detector
- DoG/LoG detector

New generation detectors:

- Stardist
- Cellpose
- Weka
- Ilastik
- MorpholibJ
- Custom segmentation via label images




TrackMate

2 - SPOT FILTERING

Filter out spots based on **spot features** such as size, average intensity, ...


Possibility to **visualize spot features** with heat maps.



TrackMate

3 - DATA ASSOCIATION

Trackers:




- **Manual tracking**
- **Overlap tracker:** based on intersection over union between two consecutive frames – only 2D
- **Nearest-neighbor tracker:** associate closest particles between two consecutive frames
- **Kalman tracker:** based on Kalman filtering, suited for objects with directed and constant motion, allows gap closing
- **Lap tracker:** implementation of utrack (Jaqaman *et al.*), based on the Minear Assignment Problem mathematical framework, allows to change linking cost between particles by penalizing features and to perform splitting, merging and gap closing
- **Simple LAP tracker:** same as LAP tracker, without cost penalties, splitting and merging

TrackMate

4 - TRACK FILTERING

Filter out tracks based on **track features** such as number of spots, duration, ...


Possibility to **visualize track features** with heat maps.



Hands on TrackMate


Open **CelegansEarly_MIP**
with **TrackMate** and
segment cells with
Thresholding detector

Filter out as many as
possible **non-cell objects**
but keep **all cells**
detected




Hands on TrackMate

Use **Simple LAP tracker**,
filter out track(s)
corresponding to non-cell
objects and **manually**
split the tracks when **cells**
divide with TrackScheme



Hands on TrackMate

Plot features, save results...



Hands on TrackMate

Use LAP tracker with
splitting and **filter out**
track(s) to directly obtain
the **2 tracks** for cells

